

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

2D HPE(TM) Large Gel NF 12.5% Kit	Kat.Nr./Cat.No. : 43304
	Lot/Contr.No. : 151377-150966

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält 4 an Kunststoffolie gebunde Gele, einschließlich Lauf- und Äquilibrationpuffer, FS Wicks und Kühlkontaktflüssigkeit. Format: 255 x 200 mm, Dicke: 0.65 mm		
	Kit comprising 4 plastic-backed gels, including running and equilibration buffers, FS wicks and cooling contact fluid. Size: 255 x 200 x 0.65 mm		
Anwendungstest application test	SDS-Elektrophorese SDS-Electrophoresis	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			04.05.2017
Lagerung (°C) storage	Buffer Kit for 2D HPE(TM) Gels		+15 bis/ to +30
Lagerung (°C) storage	2D HPE(TM) Large Gel NF-12.5 %		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH	Druckdatum/ printing date: 21.12.2015
Qualitätskontrolle/ Quality Control	

Dipl.-Ing. (FH) Bernhard Göckel	Daniela Lux-Helmstetter
---------------------------------	-------------------------

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.