

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

|                                          |                         |                      |
|------------------------------------------|-------------------------|----------------------|
| <b>2D HPE(TM) Large Gel NF 12.5% Kit</b> | <b>Kat.Nr./Cat.No.:</b> | <b>43304</b>         |
|                                          | <b>Lot/Contr.No.:</b>   | <b>160427-160324</b> |

| <b>Parameter</b><br>parameter                   | <b>Methode</b><br>method                                                                                                                                                                                                                                                                                                            | <b>Spezifikation</b><br>specification | <b>Ergebnis</b><br>result |
|-------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|---------------------------|
| <b>Beschreibung</b><br><br>description          | Kit enthält 4 an Kunststoffolie gebunde Gele, einschließlich Lauf- und Äquilibrationspuffer, FS Wicks und Kühlkontaktfüssigkeit.<br>Format: 255 x 200 mm, Dicke: 0.65 mm<br><br>Kit comprising 4 plastic-backed gels, including running and equilibration buffers, FS wicks and cooling contact fluid.<br>Size: 255 x 200 x 0.65 mm |                                       |                           |
| <b>Anwendungstest</b><br>application test       | SDS-Elektrophorese<br>SDS-Electrophoresis                                                                                                                                                                                                                                                                                           | entspricht<br>corresponds             | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life |                                                                                                                                                                                                                                                                                                                                     |                                       | 28.10.2017                |
| <b>Lagerung (°C)</b><br>storage                 | Buffer Kit for 2D HPE(TM) Gels                                                                                                                                                                                                                                                                                                      |                                       | +15 bis/ to +30           |
| <b>Lagerung (°C)</b><br>storage                 | 2D HPE(TM) Large Gel NF-12.5 %                                                                                                                                                                                                                                                                                                      |                                       | +2 bis/ to +8             |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 04.05.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.