

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Buffer Kit for 2D HPE(TM) Gels	Kat.Nr./Cat.No.:	43312
	Lot/Contr.No.:	161170

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Kit enthält Anoden- und Kathodenpuffer, Äquilibrationspuffer, Kühlkontaktflüssigkeit und Elektrodendochte. Kit comprising anode, cathode and equilibration buffer, cooling contact fluid and electrode wicks.		
Mindesthaltbarkeit minimum shelf life			04.07.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.12.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.