

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

IgG Sample Diluter IEF	Kat.Nr./Cat.No. :	43336
	Lot/Contr.No. :	130585

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare rote Lösung clear red solution	entspricht corresponds
Anwendung application	Verdünnungspuffer für CSF-Analyse in IEF-Gelen Dilution buffer required for CSF analysis in IEF gels	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		03/ 2015
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.05.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.