

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

IgG Sample Diluter IEF

Kat.Nr./Cat.No. : 43336

Lot/Contr.No. : 221183

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare rote Lösung clear red solution	entspricht corresponds
Anwendung application	Verdünnungspuffer für CSF-Analyse in IEF-Gelen dilution buffer required for CSF analysis in IEF gels	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		29.11.2024
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.12.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.