

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

IgG Sample Diluter IEF	Kat.Nr./Cat.No. :	43336
	Lot/Contr.No. :	230403

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare rote Lösung clear red solution	entspricht corresponds
Anwendung application	Verdünnungspuffer für CSF-Analyse in IEF-Gelen dilution buffer required for CSF analysis in IEF gels	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		27.03.2025
Lagerung (°C) storage		+2 bis/ to +8

SERVA Electrophoresis GmbH

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.04.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.