

# ANALYSENZERTIFIKAT

## CERTIFICATE OF ANALYSIS

**IgG Sample Diluter IEF**

**Kat.Nr./Cat.No. : 43336**

**Lot/Contr.No. : 240173**

| Parameter<br>parameter                          | Spezifikation<br>specification                                                                          | Ergebnis<br>result        |
|-------------------------------------------------|---------------------------------------------------------------------------------------------------------|---------------------------|
| <b>Aussehen</b><br>appearance                   | klare rote Lösung<br>clear red solution                                                                 | entspricht<br>corresponds |
| <b>Anwendung</b><br>application                 | Verdünnungspuffer für CSF-Analyse in IEF-Gelen<br>dilution buffer required for CSF analysis in IEF gels | entspricht<br>corresponds |
| <b>Mindesthaltbarkeit</b><br>minimum shelf life |                                                                                                         | 06.02.2026                |
| <b>Lagerung (°C)</b><br>storage                 |                                                                                                         | +2 bis/ to +8             |

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 20.02.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.