

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

CleanGel IEF for PhastSystem	Kat.Nr./Cat.No. : 43350
Size: 50 mm x 42 mm x 0.43 mm	Lot/Contr.No. : 220254

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung	20 rehydratisierbare Polyacrylamid-Gele kovalent an eine Kunststoffolie gebunden, 5 % T, ohne Probetaschen. Format: 50 x 42 mm, Dicke: 0.43 mm		
description	20 rehydratable polyacrylamide-gels covalently bound to a plastic sheet, 5 % T, without slots. Format: 50 x 42 mm, Thickness: 0.43 mm		
Mindesthaltbarkeit			11.11.2023
minimum shelf life			
Lagerung (°C)			-15 bis/ to -25
storage			

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.03.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.