

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>SERVA IEF Gel 3-9 for PhastSystem</b>	<b>Kat.-Nr./Cat. no.:</b>	<b>43366</b>
Size: 50 mm x 42 mm	<b>Lot/Contr.No. :</b>	<b>220679</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Isoelektrische Fokussierung von SERVA Kat.-Nr. 39212 Isoelectric focusing of SERVA cat. no. 39212	entspricht corresponds
<b>Inhalt</b> content	10 gebrauchsfertige IEF Gele 3-9 für das PhastSystem (Format: 50 x 42 mm) 10 ready-to use IEF gels 3-9 for the PhastSystem (size: 50 x 42 mm)	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life		04.07.2023
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 19.07.2022**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.