

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

SERVA IEF Gel 3-9 for PhastSystem	Kat.-Nr./Cat. no.:	43366
Size: 50 mm x 42 mm	Lot/Contr.No. :	220873

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Isoelektrische Fokussierung von SERVA Kat.-Nr. 39212 Isoelectric focusing of SERVA cat. no. 39212	entspricht corresponds
Inhalt content	10 gebrauchsfertige IEF Gele 3-9 für das PhastSystem (Format: 50 x 42 mm) 10 ready-to use IEF gels 3-9 for the PhastSystem (size: 50 x 42 mm)	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		23.09.2023
Lagerung (°C) storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH **Druckdatum/ printing date: 12.10.2022**
Qualitätskontrolle/ Quality Control

Dr. Judith Koch Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.