

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---------------------|--------------------------------|
| SERVA Purple | Kat.Nr./Cat.No. : 43386 |
| | Lot/Contr.No. : P130317 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|---|---|---------------------------|
| Aussehen appearance | | klare, gelbe Lösung clear, yellow solution | entspricht corresponds |
| Proteindetektion protein detection | 10x in H ₂ O mit Albumin 10x in H ₂ O with albumin | entspricht corresponds | entspricht corresponds |
| Mindesthaltbarkeit minimum shelf life | | | 07/ 2015 |
| Lagerung (°C) storage | | | -15 bis/ to -25 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 01.08.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.