

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

CSF Analysis Kit for PhastSystem	Kat.Nr./Cat.No. :	43393
	Lot/Contr.No. :	180556

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Auftrennung von CSF-Proben separation of CSF samples	entspricht corresponds
Inhalt content	CleanGel IEF for PhastSystem IEF Mix for CSF Analysis IEF sample application strips 8 slots	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		10.10.2019
Lagerung (°C) storage		-15 bis/ to -25

SERVA Electrophoresis GmbH

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.05.2018

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.