

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS**SERVA HPE™ Coomassie® Staining Kit****Kat.Nr./Cat.No. : 43396****Lot/Contr.No. : P110445**

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Sensitive Kolloidalfärbung von 1D- und 2D-Gelen nach der Elektrophorese. Highly sensitive colloidal staining of 1D and 2D gels after electrophoresis.	entspricht corresponds
Inhalt content	1000 ml Solution A 25 ml Solution B	
Mindesthaltbarkeit minimum shelf life		10/ 2013
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control**Druckdatum/ printing date: 15.08.2012**

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.