

**ANALYSENZERTIFIKAT**  
 CERTIFICATE OF ANALYSIS

<b>SERVA HPE™ Coomassie® Staining Kit</b>	<b>Kat.Nr./Cat.No. : 43396</b>
	<b>Lot/Contr.No. : P160283</b>

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	Sensitive Kolloidalfärbung von 1D- und 2D-Gelen nach der Elektrophorese. Highly sensitive colloidal staining of 1D and 2D gels after electrophoresis.	entspricht corresponds
<b>Inhalt</b> content	1000 ml Solution A 25 ml Solution B	
<b>Mindesthaltbarkeit</b> minimum shelf life		01.04.2019
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
 We do not guarantee that the product can be used for a special application.  
 This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 05.06.2018**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.