

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS



SERVA CSF Silver Staining Kit

Kat.Nr./Cat.No. : 43398

Lot/Contr.No. : P230369

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Hochsensitiver Silberfärbekit für das Färben von bis zu 5 foliengestützten, ultra-dünnen Horizontal-IEF-Gelen für die CSF-Analyse. High sensitive silver staining kit for staining of up to 5 film-backed, ultra-thin horizontal IEF gels for CSF analysis.	entspricht corresponds
Inhalt content	Fixing Solution I Ethanol denatured for Fixing Solution II and Wash Solution Fixing Solution III Solution A Solution B Solution C Solution D Solution E	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		31.03.2025
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.10.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.