

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

SERVA CSF Silver Staining Kit	Kat.Nr./Cat.No. : 43398
	Lot/Contr.No. : P240011

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	Hochsensitiver Silberfärbekit für das Färben von bis zu 5 foliengestützten, ultra-dünnen Horizontal-IEF-Gelen für die CSF-Analyse. High sensitive silver staining kit for staining of up to 5 film-backed, ultra-thin horizontal IEF gels for CSF analysis.	entspricht corresponds
Inhalt content	Fixing Solution I Ethanol denatured for Fixing Solution II and Wash Solution Fixing Solution III Solution A Solution B Solution C Solution D Solution E	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		31.03.2025
Lagerung (°C) storage		+15 bis/ to +30

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 10.01.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.