

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA PRiME™ Lightning Red	Kat.Nr./Cat.No. :	43402
	Lot/Contr.No. :	180688

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Komponenten components	SERVA PRiME™ Lightning Red DMSO, wasserfrei/ anhydrous Markierungspuffer/ Labeling Buffer Probenpuffer/ Loading Buffer		
Anwendungstest application test	SDS-PAGE	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			26.03.2020
Lagerung (°C) storage			+2 bis / to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.07.2018

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.