

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

SERVA PRiME™ Lightning Red	Kat.Nr./Cat.No. :	43402
	Lot/Contr.No. :	190077

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Komponenten components	SERVA PRiME™ Lightning Red DMSO, wasserfrei/ anhydrous Markierungspuffer/ Labeling Buffer Probenpuffer/ Loading Buffer		
Anwendungstest application test	SDS-PAGE	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			26.03.2020
Lagerung (°C) storage			+2 bis / to +8

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.01.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.