

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

SERVA PRiME™ Lightning Red	Kat.Nr./Cat.No. :	43402
	Lot/Contr.No. :	200423

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Komponenten components	SERVA PRiME™ Lightning Red DMSO, wasserfrei/ anhydrous Markierungspuffer/ Labeling Buffer Probenpuffer/ Loading Buffer		
Anwendungstest application test	SDS-PAGE	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			25.03.2022
Lagerung (°C) storage			+2 bis / to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 24.03.2021

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.