

**A N A L Y S E N Z E R T I F I K A T**  
CERTIFICATE OF ANALYSIS

---

<b>SERVA PRiME™ Lightning Red</b>	<b>Kat.Nr./Cat.No. :</b>	<b>43402</b>
	<b>Lot/Contr.No. :</b>	<b>240356</b>

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Komponenten</b> components	SERVA PRiME™ Lightning Red DMSO, wasserfrei/ anhydrous Markierungspuffer/ Labeling Buffer Probenpuffer/ Loading Buffer		
<b>Anwendungstest</b> application test	SDS-PAGE	entspricht corresponds	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			31.12.2024
<b>Lagerung (°C)</b> storage			+2 bis / to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 17.04.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.