

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Blank FocusGel 24S	Kat.Nr./Cat.No.:	43413
Size: 250 x 115 x 0.65 mm	Lot/Contr.No.:	171032

Parameter	Methode	Spezifikation	Ergebnis
parameter	method	specification	result

Beschreibung	Dieses Produkt enthält keinen SERVALYT und muss vor Gebrauch entsprechend der Produktinformation, mit SERVALYT äquilibriert werden. Die Polyacrylamid-Gele sind an eine Kunststoffolie gebunden. Format: 250 x 115 mm, Dicke: 0.65 mm, 24 Probetaschen
---------------------	---

description	This product contains no SERVALYT and has to be equilibrated with SERVALYT before use according to the product information. The polyacrylamide gels are bound to a plastic film. Format: 250 x 115 mm, Thickness: 0.65 mm, 24 sample wells
-------------	---

Mindesthaltbarkeit	07.11.2019
minimum shelf life	

Lagerung (°C)	+2 bis/ to +8
storage	

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 22.11.2017

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.