

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



Blank FocusGel 24S	Kat.Nr./Cat.No. :	43413
Size: 250 x 115 x 0.65 mm	Lot/Contr.No. :	230834

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung description	Dieses Produkt enthält keinen SERVALYT und muss vor Gebrauch entsprechend der Produktinformation, mit SERVALYT rehydratisiert werden. Die Polyacrylamid-Gele sind an eine Kunststoffolie gebunden. Format: 250 x 115 mm, Dicke: 0.65 mm, 24 Probetaschen This product contains no SERVALYT and has to be rehydrated with SERVALYT before use according to the product information. The polyacrylamide gels are bound to a plastic film. Format: 250 x 115 mm, Thickness: 0.65 mm, 24 sample wells		
Mindesthaltbarkeit minimum shelf life	20.07.2025		
Lagerung (°C) storage	+2 bis/ to +8		

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 03.08.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.