

ANALYSENZERTIFIKAT
 CERTIFICATE OF ANALYSIS

SERVA SDS Gel 8-25 % Kit for PhastSystem	Kat.-Nr./Cat. No.:	43503
Size: 50 mm x 42 mm	Lot/Contr.No. :	230125

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung	10 gebrauchsfertige SDS-Polyacrylamid-Gele kovalent an eine Kunststoffolie gebunden (Format: 50 x 42 mm) 20 Elektroden-Wicks SDS-Kathoden- und SDS-Anoden-Puffer (je 20 ml)		
description	10 ready-to-use SDS polyacrylamide gels covalently bound to a plastic sheet (size: 50 x 42 mm) 20 electrode wicks SDS cathode and SDS anode buffer (20 ml each)		
Mindesthaltbarkeit minimum shelf life			24.01.2024
Lagerung (°C) storage			+ 2 bis/ to + 8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH **Druckdatum/ printing date: 01.02.2023**
Qualitätskontrolle/ Quality Control

Dr. Judith Koch Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.