

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



| | | |
|---|---------------------------|---------------|
| SERVA SDS Gel 8-25 % Kit for PhastSystem | Kat.-Nr./Cat. No.: | 43503 |
| Size: 50 mm x 42 mm | Lot/Contr.No. : | 231230 |

| Parameter parameter | Methode method | Spezifikation specification | Ergebnis result |
|---|---|---|---------------------------|
| Beschreibung description | 10 gebrauchsfertige SDS-Polyacrylamid-Gele kovalent an eine Kunststoffolie gebunden (Format: 50 x 42 mm) 20 Elektroden-Wicks SDS-Kathoden- und SDS-Anoden-Puffer (je 20 ml) | 10 ready-to-use SDS polyacrylamide gels covalently bound to a plastic sheet (size: 50 x 42 mm) 20 electrode wicks SDS cathode and SDS anode buffer (20 ml each) | |
| Mindesthaltbarkeit minimum shelf life | | | 04.12.2025 |
| Lagerung (°C) storage | | | + 2 bis/ to + 8 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.12.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.