

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA SDS Gel 8-25 % Kit for PhastSystem	Kat.-Nr./Cat. No.:	43503
Size: 50 mm x 42 mm	Lot/Contr.No. :	240884

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Beschreibung	10 gebrauchsfertige SDS-Polyacrylamid-Gele kovalent an eine Kunststoffolie gebunden (Format: 50 x 42 mm) 20 Elektroden-Wicks SDS-Kathoden- und SDS-Anoden-Puffer (je 20 ml)		
description	10 ready-to-use SDS polyacrylamide gels covalently bound to a plastic sheet (size: 50 x 42 mm) 20 electrode wicks SDS cathode and SDS anode buffer (20 ml each)		
Mindesthaltbarkeit minimum shelf life	27.09.2025		
Lagerung (°C) storage	+ 2 bis/ to + 8		

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 02.10.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.