

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS



SERVA SDS Gel 8-25 % Kit for PhastSystem	Kat.-Nr./Cat. No.:	43503
Size: 50 mm x 42 mm	Lot/Contr.No. :	250094

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
-------------------------------	--------------------------	---------------------------------------	---------------------------

Beschreibung
10 gebrauchsfertige SDS-Polyacrylamid-Gele kovalent an eine Kunststoffolie gebunden (Format: 50 x 42 mm)
20 Elektroden-Wicks
SDS-Kathoden- und SDS-Anoden-Puffer (je 20 ml)

description
10 ready-to-use SDS polyacrylamide gels covalently bound to a plastic sheet (size: 50 x 42 mm)
20 electrode wicks
SDS cathode and SDS anode buffer (20 ml each)

Mindesthaltbarkeit
minimum shelf life
27.09.2025

Lagerung (°C)
storage
+ 2 bis/ to + 8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.01.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.