

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Visking® Dialysierschlauch	Kat.Nr./Cat.No. :	44120
Visking® dialysis tubing	Lot/Contr.No. :	251081

Parameter parameter	Spezifikation specification	Ergebnis result
Material	regenerierte Cellulose regenerated cellulose	entspricht corresponds
Flachbreite (mm) flat width	44	entspricht corresponds
Durchmesser (mm) diameter	28	entspricht corresponds
Porendurchmesser (Å) pore diameter	ca. 25	entspricht corresponds
Ausschlußvolumen (Dalton) exclusion limit	12 000 – 14 000	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		29.10.2030
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 29.10.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.