

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Visking® Dialysierschlauch	Kat.Nr./Cat.No. : 44130
Visking® dialysis tubing	Lot/Contr.No. : 170941

Parameter parameter	Spezifikation specification	Ergebnis result
Material	regenerierte Cellulose regenerated cellulose	entspricht corresponds
Flachbreite (mm) flat width	117	entspricht corresponds
Durchmesser (mm) diameter	75	entspricht corresponds
Porendurchmesser (Å) pore diameter	ca. 25	entspricht corresponds
Ausschlußvolumen (Dalton) exclusion limit	12 000 – 14 000	entspricht corresponds
Mindesthaltbarkeit minimum shelf life		17.10.2022
Lagerung (°C) storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.