

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Acetonitril für HPLC	Kat.Nr./Cat.No. :	45605
Acetonitrile for HPLC	Lot/Contr.No. :	120125

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		CH ₃ CN	
Molekulargewicht molecular weight		41.05	
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Gehalt (%) assay	GC	min. 99.8	99.99
Brechungsindex refractive index		1.3430 – 1.3450	1.3439
Wasser (%) water	KF	max. 0.03	0.01
T 1cm (%)	gegen HPLC-Wasser against HPLC-water	210 nm 240 nm	min. 93 min. 99 97 100
Mindesthaltbarkeit minimum shelf life			08/ 2014
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 01.02.2012

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.