

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Isopropanol p.A.	Kat.Nr./Cat.No. :	45629
Isopropanol	Lot/ Contr.No. :	131235

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		C ₃ H ₈ O	
Molekulargewicht molecular weight		60.1	
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Gehalt (%) assay	GC	min. 99.7	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (g/ml) specific density	20 °C	0.785	entspricht corresponds
Brechungsindex refractive index		1.3770 – 1.3780	1.3770
Wasser (ppm) water		max. 1000	entspricht corresponds
Acidität (%) acidity		max. 0.001	entspricht corresponds
Abdamprückstand (ppm) Residue on evaporation		max. 10	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			10/ 2016
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.11.2013

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.