

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Methanol p.A.	Kat.Nr./Cat.No. :	45631
Methanol analytical grade	Lot/Contr.No. :	160306

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Summenformel formula		CH ₃ OH	
Molekulargewicht molecular weight		32.04	
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Gehalt (%) assay	GC	min. 99.8	≥ 99.8
Brechungsindex refractive index	n _D ²⁰	1.328 – 1.329	1.3285
Dichte (g/ml) specific density	20 °C	0.7915	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Wasser (ppm) water	KF	max. 1000	≤ 500
Freie Säure (ppm) acidity	als/ as HCOOH	max. 20	≤ 20
Abdampfrückstand (ppm) residue on evaporation		max. 10	≤ 10
Mindesthaltbarkeit minimum shelf life			19.03.2019
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 30.03.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.