

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Colchicin Lösung sterilfiltriert	Kat.Nr./Cat.No. :	47252
Colchicin solution sterile filtered	Lot/Contr.No. :	P160041

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare farblose Lösung clear colorless solution	entspricht corresponds
Konzentration (µg/ml) concentration	in PBS	10	10
pH		7.0 – 8.0	7.3
Osmolarität (mOsmol/l) osmolarity		270 – 300	291
Sterilität sterility		negativ negative	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			29.01.2018
Lagerung (°C) storage			-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.02.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.