

CERTIFICATE OF ANALYSIS

Collagen R Solution 0.2%
sterile

Cat.No.: 47254
Contr.No. : 090595

Parameter	Method	Specification	Result
Molecular weight		ca. 300 000	
Source		rat tail	
Concentration		2 mg/ ml in 0.1% acetic acid	
Sterility		passes test	corresponds
A 1cm/ 280 nm	undiluted	max. 0.2	0.066
Gelling test		corresponds to reference	corresponds to reference
Minimum shelf life			06/ 2011
Storage (°C)			+2 to +8
References		Strom, S. C. & Michalopoulos, G. (1982) Methods Enzymol. 82 , 544 - 55 Miller, E. J. (1976) Mol. Cell. Biochem. 13 , 165 - 91	

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.04.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.