

**A N A L Y S E N Z E R T I F I K A T**  
 CERTIFICATE OF ANALYSIS

<b>G 418 Lösung 50 mg/ ml, sterilfiltriert</b>	<b>Kat.Nr./Cat.No. : 47995</b>
<b>G 418 solution 50 mg/ ml, steril filtered</b>	<b>Lot/Contr.No. : P200088</b>

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Aussehen</b> appearance		klare, farblose bis gelbliche Lösung clear, colorless to pale yellow solution	entspricht corresponds
<b>Konzentration (mg/ ml in H<sub>2</sub>O)</b> concentration		50 mg G-418	entspricht corresponds
<b>pH</b>		4.0 – 6.0	5.1
<b>Osmolarität (mOsmol/ l)</b> osmolarity			155
<b>Sterilttest</b> sterile test		kein Wachstum no growth	entspricht corresponds
<b>Mindesthaltbarkeit</b> minimum shelf life			03.03.2022
<b>Lagerung (°C)</b> storage			-15 bis/ to -25

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 09.03.2020**

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.