

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

LB, Pulvermedium	Kat.Nr./Cat.No. : 48501
LB medium, powder	Lot/Contr.No. : P240387

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		elfenbeinfarbenes Pulver ivory-colored powder	entspricht corresponds
Löslichkeit solubility	20 g/l in H ₂ O	goldgelbe Lösung gold-yellowish solution	entspricht corresponds
pH	20 g/l in H ₂ O	7.0 ± 0.3	6.8
Osmolarität (mOsmol/l) osmolarity	20 g/l in H ₂ O		239
Zusammensetzung composition			
Tryptone (# 48647)		10.0 g/l	entspricht corresponds
Yeast extract (# 24540)		5.0 g/l	entspricht corresponds
NaCl (# 30183)		5.0 g/l	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			02.12.2027
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.12.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.