

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

2x YT Pulvermedium	Kat.Nr./Cat.No. : 48503
2x YT medium, powder	Lot/Contr.No. : P140161

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes Pulver light ivory-colored powder	entspricht corresponds
Löslichkeit solubility	31 g/l in H ₂ O	klare gelbe Lösung clear yellow solution	entspricht corresponds
pH	31 g/l in H ₂ O	6.5 – 7.5	7.0
Osmolarität (mOsmol/l) osmolarity	31 g/l in H ₂ O		291
Zusammensetzung composition			
Tryptone (# 48647)		16.0 g/l	entspricht corresponds
Yeast extract (# 24540)		10.0 g/l	entspricht corresponds
NaCl (# 30183)		5.0 g/l	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			26.03.2018
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 04.05.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.