

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

2x YT Pulvermedium	Kat.Nr./Cat.No. : 48503
2x YT medium, powder	Lot/Contr.No. : P190045

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellelfenbeinfarbenes Pulver light ivory-colored powder	entspricht corresponds
Löslichkeit solubility	31 g/l in H ₂ O	klare gelbe Lösung clear yellow solution	entspricht corresponds
pH	31 g/l in H ₂ O	6.5 – 7.5	6.9
Osmolarität (mOsmol/l) osmolarity	31 g/l in H ₂ O		289
Zusammensetzung composition			
Tryptone (# 48647)		16.0 g/l	entspricht corresponds
Yeast extract (# 24540)		10.0 g/l	entspricht corresponds
NaCl (# 30183)		5.0 g/l	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			23.01.2022
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
 Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.01.2019

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.