

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Sojabohne	Kat.Nr./Cat.No. :	57556
Lecithin from soybean	Lot/Contr.No. :	130571

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellbeiges Pulver light beige powder	entspricht corresponds
Acetonunlös. Substanz (%) acetone insol. substance		min. 95.0	96.9
Toluollösliche Substanzen (%) toluol soluble substances		min. 97.5	≥ 99.7
Säurezahl (mg KOH/ g) acid number		max. 30.0	23.7
Phosphatidylcholin (%) phosphatidylcholine		18.0 – 25.0	20.7
Phosphatidylethanolamin (%) phosphatidylethanolamine		10.0 – 16.0	12.4
Phosphatidylinositol (%)		14.0 – 19.0	16.0
Phosphatidsäure (%) phosphatic acid		4.0 – 7.0	5.0
pH-Wert pH-value	1 % in H ₂ O	6.0 – 7.0	6.9
Wasser (%) water		max. 1.5	1.23
Mindesthaltbarkeit minimum shelf life			10/ 2014
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 16.05.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.