

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Sojabohne	Kat.Nr./Cat.No. :	57556
Lecithin from soybean	Lot/Contr.No. :	140493

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellbeiges Pulver light beige powder	entspricht corresponds
Acetonunlös. Substanz (%) acetone insol. substance		min. 96.5	96.8
Toluolunlösliche Substanzen (%) toluol insoluble substances		max. 0.3	entspricht corresponds
Säurezahl (mg KOH/ g) acid number		max. 30.0	25.5
Phosphatidylcholin (%) phosphatidylcholine		17.0 – 25.0	18.3
Phosphatidylethanolamin (%) phosphatidylethanolamine		10.0 – 16.0	11.1
Phosphatidylinositol (%)		13.0 – 19.0	14.4
Phosphatidsäure (%) phosphatic acid		2.0 – 9.0	6.0
pH-Wert pH-value	1 % in H ₂ O	6.0 – 7.0	6.7
Wasser (%) water		max. 1.5	1.05
Mindesthaltbarkeit minimum shelf life			12/ 2015
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.04.2014

Dipl.-Ing. (FH) Bernhard Göckel

Jakob Arndt

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.