

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Lecithin aus Sojabohne	Kat.Nr./Cat.No. :	57556
Lecithin from soybean	Lot/Contr.No. :	140977

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellbeiges Pulver lightbeige powder	entspricht corresponds
Acetonunlös. Substanz (%) acetone insol. substance		min. 96.5	96.8
Toluolunlösliche Substanzen (%) toluol insoluble substances		max. 0.3	entspricht corresponds
Säurezahl (mg KOH/ g) acid number		max. 30.0	25.5
Phosphatidylcholin (%) phosphatidylcholine		17.0 – 25.0	18.3
Phosphatidylethanolamin (%) phosphatidylethanolamine		10.0 – 16.0	11.1
Phosphatidylinositol (%)		13.0 – 19.0	14.4
Phosphatidsäure (%) phosphatic acid		2.0 – 9.0	6.0
pH-Wert pH-value	1 % in H ₂ O	6.0 – 7.0	6.7
Wasser (%) water		max. 1.5	1.1
Mindesthaltbarkeit minimum shelf life			31.12.2016
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.10.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.