

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

Lecithin aus Sojabohne	Kat.Nr./Cat.No. :	57556
Lecithin from soybean	Lot/Contr.No. :	160693

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellbeiges Pulver lightbeige powder	entspricht corresponds
Acetonunlös. Substanz (%) acetone insol. substance		min. 96.5	96.8
Toluolunlösliche Substanzen (%) toluol insoluble substances		max. 0.3	entspricht corresponds
Säurezahl (mg KOH/ g) acid number		max. 30.0	26.9
Phosphatidylcholin (%) phosphatidylcholine		17.0 – 25.0	18.9
Phosphatidylethanolamin (%) phosphatidylethanolamine		10.0 – 16.0	12.5
Phosphatidylinositol (%)		13.0 – 19.0	15.1
Phosphatidsäure (%) phosphatic acid		2.0 – 9.0	5.9
pH-Wert pH-value	1 % in H ₂ O	6.0 – 7.0	6.6
Wasser (%) water		max. 1.5	0.9
Mindesthaltbarkeit minimum shelf life			31.01.2018
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 19.07.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.