

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

Lecithin aus Sojabohne	Kat.Nr./Cat.No. :	57556
Lecithin from soybean	Lot/Contr.No. :	170292

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		hellbeiges Pulver lightbeige powder	entspricht corresponds
Acetonunlös. Substanz (%) acetone insol. substance		min. 96.5	97.1
Toluolunlösliche Substanzen (%) toluol insoluble substances		max. 0.3	entspricht corresponds
Säurezahl (mg KOH/ g) acid number		max. 30.0	25.1
Phosphatidylcholin (%) phosphatidylcholine		17.0 – 25.0	19.7
Phosphatidylethanolamin (%) phosphatidylethanolamine		10.0 – 16.0	12.5
Phosphatidylinositol (%)		13.0 – 19.0	15.7
Phosphatidsäure (%) phosphatic acid		2.0 – 9.0	4.6
pH-Wert pH-value	1 % in H ₂ O	6.0 – 7.0	6.9
Wasser (%) water		max. 1.5	0.8
Mindesthaltbarkeit minimum shelf life			28.02.2019
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 11.04.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.